
maXXi 5

Röntgenfluoreszenz Analysator X-ray Fluorescence Analyser

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A universal X-ray coating
thickness and material
analyser designed for the
inspection of large and
small precision parts

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The maXXi 5 offers a universal coating thickness and material analyser at a reasonable price/ performance ratio. Intended for customers where in addition to highest precision ease of use and versatility are important the maXXi 5 is the ideal measuring system. Precise, non-destructive measurement is achieved in only a few seconds. It is possible to determine not only the thickness but also the composition of alloy coatings. In addition the maXXi 5 can also be used for the analysis of plating solutions.

The voluminous chamber together with optimised X-ray beam geometry allows the use of very small collimators. A multiple position collimator changer in conjunction with programmable X-ray tube power enables adaptation for almost any required application.

A software controlled X-Y-Z-sample stage allows automatic measurement sequences. This makes the maXXi 5 an ideal tool for the measurement of virtually all metallic coatings.

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

Technical Data

HV Generator	50 kV, max. 1.2 mA current, software controlled
X-ray tube	Micro focus tube, W-target, Spot 85 µm x 85 µm, glass window
X-ray power	max. 50 VA, applications optimised
(Option)	Micro focus tube, spot 85 µm x 85 µm, Be-window
Single collimator	i. e. 0.1 mm or 0.3 mm or 0.1 mm x 0.3 mm
Collimator changer	Four positions, programmable and motorised
(Option)	from 60 x 60 µm to 300 µm Ø or customs specified
Sample Stage (Option)	a) Programmable and motorized X 300 mm, Y 250 mm, Z 180 mm, Easy Load, max. load 5 kg b) Programmable and motorized X 300 mm, Y 250 mm, Z 180 mm, max. load 20 kg
Sample positioning	Joy Stick, Point & Shoot, Laser, Auto-focus
Video system	Colour video system, field of view 2 x 1.5 mm and 4 x 3 mm
Magnification	approx. 40 x and 80 x
Radiation safety	PTB tested, totally fail safe system
Chamber	H 67 cm, B 60 cm, D 75 cm, slotted sample chamber
Sample Chamber	Useable max. H 35 cm, W 60 cm, D 50 cm
Mains supply	110 V / 230 V 60 Hz / 50 Hz
X-Master-Software	Operation system with WINDOWS 2000 / XP using modern PC Technology
µ-Master	Evaluation module for coating thickness measurement
Fun-Master	Calibration module for fundamental parameter mathematics
Element-Master	Standardless qualitative and quantitative material analysis for up to 20 elements
%-Master	Precise quantitative material analysis for up to 8 elements simultaneously with statistics
Liquid-Master	Plating bath analysis
Data-Master	Data base and statistic software with long term documentation
Report-Master	Software module to create customised reports

All specifications are subject to change without notice

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GENESEEA TECHNOLOGY CO., LTD.
SHENZHEN GENESEA TECHNOLOGY CO., LTD.